

# Recombinant Hantaan virus HTNV (strain 84FLi) Nucleocapsid / NP Protein (His Tag)

Catalog Number:PKSV030173

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


## Description

<b>Synonyms</b>	NP Protein;HTNV
<b>Species</b>	HTNV
<b>Expression Host</b>	Baculovirus-Insect Cells
<b>Sequence</b>	Met1-Leu429
<b>Accession</b>	AAK66740.1
<b>Calculated Molecular Weight</b>	49.7 kDa.
<b>Tag</b>	C-His

## Properties

<b>Purity</b>	> 80 % as determined by reducing SDS-PAGE.
<b>Endotoxin</b>	< 1.0 EU per $\mu$ g of the protein as determined by the LAL method.
<b>Storage</b>	Store at $-20^{\circ}\text{C}$ , stable for 6 months. Please minimize freeze-thaw cycles.
<b>Shipping</b>	This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at $-20^{\circ}\text{C}$ .
<b>Formulation</b>	Supplied as sterile 20 mM Tris, 500 mM NaCl, 20 % glycerol, pH 8.0.
<b>Reconstitution</b>	Not Applicable

## Data



> 80 % as determined by reducing SDS-PAGE.

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017