

Recombinant WNV (lineage 1, strain NY99) NS1 Protein (His Tag)



Catalog Number:PKSV030261

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	WNV-NS1 Protein;WNV;West Nile Virus;strain NY99
Species	WNV
Expression Host	HEK293 Cells
Sequence	Asp1-Ala352
Accession	YP_001527881.1
Calculated Molecular Weight	42.1 kDa.
Tag	N-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per μg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS,pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

interferon regulatory factor 3

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017