

RBMX2 Polyclonal Antibody

Catalog Number:E-AB-66269

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Mouse
Immunogen	Recombinant fusion protein of human RBMX2 (NP_057108.2).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB	1:500-1:2000
-----------	--------------

Data



Western blot analysis of extracts of Mouse liver using RBMX2 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:37kDa
Calculated Mw:37kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

RBMX2 (RNA Binding Motif Protein X-Linked 2) is a Protein Coding gene. Diseases associated with RBMX2 include Spermatogenic Failure, X-Linked, 1. Gene Ontology (GO) annotations related to this gene include nucleic acid binding and nucleotide binding. An important paralog of this gene is CSTF2T.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017