

CFAP61 Polyclonal Antibody

Catalog Number:E-AB-67314

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse
Immunogen	A synthetic peptide of human CFAP61.
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB	1:500-1:2000
-----------	--------------

Data



Western blot analysis of extracts of Mouse eye using CFAP61 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:141KDa
Calculated Mw:141kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

CFAP61 (Cilia And Flagella Associated Protein 61) is a Protein Coding gene. Diseases associated with CFAP61 include Rigid Spine Muscular Dystrophy 1 and Corneal Dystrophy, Endothelial, X-Linked. Gene Ontology (GO) annotations related to this gene include oxidoreductase activity.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017