

STMN4 Polyclonal Antibody

Catalog Number:E-AB-67534



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Mouse,Rat
Immunogen	Recombinant fusion protein of human STMN4 (NP_110422.2).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB	1:500-1:1000
-----------	--------------

Data



Western blot analysis of extracts of various cell lines using STMN4 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:33kDa

Calculated Mw:19kDa/22kDa/23kDa/25kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

STMN4 (Stathmin 4) is a Protein Coding gene. Diseases associated with STMN4 include Spastic Paraplegia 36, Autosomal Dominant. Gene Ontology (GO) annotations related to this gene include tubulin binding. An important paralog of this gene is STMN3.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017