

PPP3CB Polyclonal Antibody

Catalog Number:E-AB-67584



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse
Immunogen	Recombinant fusion protein of human PPP3CB (NP_001135825.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB	1:500-1:2000
-----------	--------------

Data



Western blot analysis of extracts of Mouse brain using PPP3CB Polyclonal Antibody at dilution of 1:1000.

Observed Mw:69kDa

Calculated Mw:58kDa/59kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

PPP3CB (Protein Phosphatase 3 Catalytic Subunit Beta) is a Protein Coding gene. Diseases associated with PPP3CB include Pancreatic Ductal Adenocarcinoma and Kidney Hypertrophy. Among its related pathways are Nur77 Signaling in T-Cell and B cell receptor signaling pathway (KEGG). Gene Ontology (GO) annotations related to this gene include calcium ion binding and hydrolase activity. An important paralog of this gene is PPP3CA.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017