

Monoclonal Mouse Anti-Human IgE(C3-domain specific)(Biotin conjugated)

Catalog No. E-AB-1073

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Human
Immunogen	Human IgE C3 domain
Host	Mouse
Clone	4F4cc
Purification	Affinity purification
Conjugation	Biotin
Buffer	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

Applications Recommended Dilution

Western Blotting	1:200-4000(Enhanced chemiluminescent detection)
Western Blotting	1:200-400(DAB detection)
Immunohistochemistry	1:100-200
Direct ELISA	1:1000-5000(TMB detection)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.