

Monoclonal Mouse Anti-Human IgE(C3-domain specific)(Biotin conjugated)



Catalog Number:E-AB-1073

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Human
Immunogen	Human IgE C3 domain
Host	Mouse
Clone	4F4cc
Purification	Affinity purification
Conjugation	Biotin
Formulation	PBS with 1% protective protein, 0.04% Proclin 300 and 50% glycerol, PH7.4

Applications Recommended Dilution

Western Blotting	1:200-4000(Enhanced chemiluminescent detection)
Western Blotting	1:200-400(DAB detection)
Immunohistochemistry	1:100-200
Direct ELISA	1:1000-5000(TMB detection)

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017