

Acetyl Lysine Monoclonal Antibody

Catalog No. E-AB-20301

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	All
Immunogen	Purified Protein
Host	Mouse
Isotype	IgG
Clone	Clone:8E3
Purification	Protein A purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.


Applications Recommended Dilution

WB	1:1000-2000
IHC	1:50-1:200

Data



Western Blot analysis of 1) Mouse brain, 2) HeLa, 3) HeLa+TSA treated with Acetyl Lysine Monoclonal Antibody.



Immunohistochemistry of paraffin-embedded Human breast carcinoma tissue with Acetyl Lysine Monoclonal Antibody.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only