

Rubisco (Large Chain) Monoclonal Antibody

Catalog No. E-AB-22122

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Plants
Immunogen	Protein
Host	Mouse
Isotype	IgG
Purification	Protein A purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 0.5% protective protein and 50% glycerol, pH7.4

Applications Recommended Dilution

WB 1:2000-5000

Data



Western Blot analysis of Arabidopsis using Rubisco (Large Chain) Monoclonal Antibody at dilution of 1)1:2000 2)1:5000.

Observed Mw:53kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only