

PDGFRA Monoclonal Antibody

Catalog No. E-AB-22143

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Human,Mouse,Rat
Immunogen	Synthetic Peptide of PDGFR α
Host	Mouse
Isotype	IgG
Clone	4E7
Purification	Protein A purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 0.5% protective protein and 50% glycerol, pH7.4


Applications Recommended Dilution

IHC	1:100-200
IF	1:100-200

Data



Immunohistochemistry of paraffin-embedded Human skin tissue using PDGFRA Monoclonal Antibody at dilution of 1:200.



Immunofluorescence analysis of Mouse spleen tissue using PDGFRA Monoclonal Antibody at dilution of 1:200.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Receptor that binds both PDGFA and PDGFB and has a tyrosine-protein kinase activity.

For Research Use Only