

Bovine Serum Albumin Monoclonal Antibody

Catalog Number:E-AB-22156



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Cow
Immunogen	Purified Protein
Host	Mouse
Isotype	IgG
Purification	Protein A purification
Formulation	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB	1:1000-2000
-----------	-------------

Data



Western Blot analysis of Bovine Serum Albumin with Bovine Serum Albumin Monoclonal Antibody.

Observed Mw:66kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017