

Rubisco (Large Chain) Monoclonal Antibody

Catalog Number:E-AB-22168



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Plants
Immunogen	Protein
Host	Mouse
Isotype	IgG
Purification	Protein A purification
Formulation	PBS with 0.02% sodium azide, 0.5% protective protein and 50% glycerol, pH7.4

Applications Recommended Dilution

WB 1:2000-5000

Data



Western Blot analysis of Arabidopsis using Rubisco (Large Chain) Monoclonal Antibody at dilution of 1) 1:2000, 2) 1:5000.

Observed Mw:53kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017