

## Donkey Anti-Rabbit IgG(H+L)(peroxidase/HRP conjugated)

Catalog No. E-AB-1080

*Note:* Centrifuge before opening to ensure complete recovery of vial contents.

### Description

<b>Immunogen</b>	Rabbit IgG
<b>Host</b>	Donkey
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	HRP
<b>Buffer</b>	PBS with 50% glycerol, PH7.4

### Applications Recommended Dilution

**Western Blotting:**  
**1:2000-5000**  
(Enhanced chemiluminescent detection ) Western Blotting:  
**1:300-1500(DAB detection) Direct**  
**ELISA:**  
**1:3000-5000(TMB detection)**

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

### For Research Use Only