

Monoclonal Mouse Anti-Human IgE(C3-domain specific)

Catalog No. E-AB-1072

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Immunogen	Human IgE C3 domain
Host	Mouse
Clone	4F4cc
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

Applications Recommended Dilution

Other When the antibody is labeled, for example, with HRP, FITC, BIOTIN, etc, BSA and NaN3 cannot be contained in reagents.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.