

## ZSCAN4C Polyclonal Antibody

**Catalog No.** E-AB-63298

*Note:* Centrifuge before opening to ensure complete recovery of vial contents.


### Description

<b>Reactivity</b>	Mouse
<b>Immunogen</b>	A synthetic peptide of mouse ZSCAN4
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	Unconjugated
<b>Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

### Applications Recommended Dilution

**WB** 1:500-1:2000

### Data



Western blot analysis of extracts of Mouse kidney using ZSCAN4C Polyclonal Antibody at dilution of 1:1000.

**Observed Mw:57kDa**  
**Calculated Mw:57kDa**

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

### Background

ZSCAN4 (Zinc Finger And SCAN Domain Containing 4) is a Protein Coding gene. Diseases associated with ZSCAN4 include Facioscapulohumeral Muscular Dystrophy 1 and Facioscapulohumeral Muscular Dystrophy 2. Among its related pathways are Preimplantation Embryo. Gene Ontology (GO) annotations related to this gene include DNA-binding transcription factor activity. An important paralog of this gene is ZNF263.

### For Research Use Only