

MFAP2 Polyclonal Antibody

Catalog No. E-AB-63308

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Recombinant protein of human MFAP2
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB 1:500 - 1:2000

Data



Western blot analysis of extracts of various cell lines with MFAP2 Polyclonal Antibody

Observed Mw:31kDa

Calculated Mw:17kDa/20kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Microfibrillar-associated protein 2 is a major antigen of elastin-associated microfibrils and a candidate for involvement in the etiology of inherited connective tissue diseases. Four transcript variants encoding two different isoforms have been found for this gene.

For Research Use Only