

Vimentin Polyclonal Antibody

Catalog No. E-AB-64077

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description


Reactivity	Human, Mouse, Rat
Immunogen	Recombinant protein of human VIM
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB 1:500 - 1:1000


IHC 1:50 - 1:100

Data




Western blot analysis of extracts of various cell lines with Vimentin Polyclonal Antibody

Observed Mw:53kDa
Calculated Mw:53kDa



Immunohistochemistry of paraffin-embedded mouse brain with Vimentin Polyclonal Antibody



Immunohistochemistry of paraffin-embedded rat brain with Vimentin Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

For Research Use Only

This gene encodes a member of the intermediate filament family. Intermediate filamentents, along with microtubules and actin microfilaments, make up the cytoskeleton. The protein encoded by this gene is responsible for maintaining cell shape, integrity of the cytoplasm, and stabilizing cytoskeletal interactions. It is also involved in the immune response, and controls the transport of low-density lipoprotein (LDL)-derived cholesterol from a lysosome to the site of esterification. It functions as an organizer of a number of critical proteins involved in attachment, migration, and cell signaling. Mutations in this gene causes a dominant, pulverulent cataract.