

# Anti-Ebola virus EBOV(subtype Sudan, strain Gulu) VP40/Matrix protein VP40 Polyclonal Antibody



Catalog Number:E-AB-V1116

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

|                     |  |
|---------------------|--|
| <b>Reactivity</b>   | Ebola virus  |
| <b>Immunogen</b>    | Recombinant EBOV (subtype Sudan, strain Gulu) VP40 / Matrix protein VP40 Protein (His Tag) |
| <b>Host</b>         | Rabbit   |
| <b>Isotype</b>      | IgG  |
| <b>Purification</b> | Protein A Affinity   |
| <b>Conjugation</b>  | Unconjugated   |
| <b>Formulation</b>  | 0.2 µm filtered solution in PBS  |

## Applications Recommended Dilution

|              |               |
|--------------|---------------|
| <b>ELISA</b> | 1:1000-1:2000 |
|--------------|---------------|

## Preparation & Storage

|                |   |
|----------------|---|
| <b>Storage</b> | Store at -20°C. Avoid freeze / thaw cycles. |
|----------------|---|

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017