

AMBRA1 Polyclonal Antibody

Catalog No. E-AB-64255

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Mouse
Immunogen	Recombinant protein of human AMBRA1
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of mouse brain with
AMBRA1 Polyclonal Antibody

Observed Mw:160kDa

Calculated

Mw:84kDa/132kDa/135kDa/139kDa/142kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only