

TDRD1 Polyclonal Antibody

Catalog No. E-AB-64357

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Mouse
Immunogen	Recombinant fusion protein of human TDRD1 (NP_942090.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:1000-1:2000

Data



Western blot analysis of extracts of Mouse testis using TDRD1 Polyclonal Antibody at dilution of 1:3000.

Observed Mw:132kDa

Calculated Mw:79kDa/119kDa/132kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

This gene is similar to a mouse gene that encodes a tudor domain protein. Alternatively spliced transcript variants have been described but their full length sequences have not been determined.

For Research Use Only