

CHCHD6 Polyclonal Antibody

Catalog No. E-AB-64373

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description


Reactivity	Human
Immunogen	Recombinant fusion protein of human CHCHD6 (NP_115719.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications

Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines using CHCHD6 Polyclonal Antibody at dilution of 1:3000.

Observed Mw:26kDa
Calculated Mw:26kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

CHCHD6 (Coiled-Coil-Helix-Coiled-Coil-Helix Domain Containing 6) is a Protein Coding gene. Diseases associated with CHCHD6 include Tinea Manuum. An important paralog of this gene is CHCHD3.

For Research Use Only