

HEATR1 Polyclonal Antibody

Catalog No. E-AB-64455

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Recombinant fusion protein of human HEATR1 (NP_060542.4).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of HeLa cells using HEATR1 Polyclonal Antibody at dilution of 1:3000.

Observed Mw:242kDa
Calculated Mw:242kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

HEATR1 (HEAT Repeat Containing 1) is a Protein Coding gene. Among its related pathways are Gene Expression and rRNA processing in the nucleus and cytosol. Gene Ontology (GO) annotations related to this gene include snoRNA binding.

For Research Use Only