

HEATR1 Polyclonal Antibody

Catalog Number:E-AB-64455



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Recombinant fusion protein of human HEATR1 (NP_060542.4).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB	1:500-1:2000
-----------	--------------

Data



Western blot analysis of extracts of HeLa cells using HEATR1 Polyclonal Antibody at dilution of 1:3000.

Observed Mw:242kDa
Calculated Mw:242kDa

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

Background

HEATR1 (HEAT Repeat Containing 1) is a Protein Coding gene. Among its related pathways are Gene Expression and rRNA processing in the nucleus and cytosol. Gene Ontology (GO) annotations related to this gene include snoRNA binding.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017