

REG4 Polyclonal Antibody

Catalog Number:E-AB-64471

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Mouse,Rat
Immunogen	Recombinant fusion protein of human REG4 (NP_114433.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB	1:500-1:2000
-----------	--------------

Data



Western blot analysis of extracts of various cell lines using REG4 Polyclonal Antibody at dilution of 1:3000.

Observed Mw:18kDa

Calculated Mw:14kDa/18kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

REG4 (Regenerating Family Member 4) is a Protein Coding gene. Diseases associated with REG4 include Ulcerative Colitis and Pseudomyxoma Peritonei. Among its related pathways are Gastric cancer and Adhesion. Gene Ontology (GO) annotations related to this gene include calcium ion binding and heparin binding. An important paralog of this gene is REG3A.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017