

## ERCC6L2 Polyclonal Antibody

**Catalog No.** E-AB-64479

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


### Description

|                     |                                                      |
|---------------------|------------------------------------------------------|
| <b>Reactivity</b>   | Human, Mouse                                         |
| <b>Immunogen</b>    | Recombinant protein of human ERCC6L2                 |
| <b>Host</b>         | Rabbit                                               |
| <b>Isotype</b>      | IgG                                                  |
| <b>Purification</b> | Affinity purification                                |
| <b>Conjugation</b>  | Unconjugated                                         |
| <b>Buffer</b>       | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. |

### Applications Recommended Dilution

**WB** 1:500-1:2000

### Data



Western blot analysis of extracts of various cell lines  
with ERCC6L2 Polyclonal Antibody

**Observed Mw:81kDa**

**Calculated Mw:60kDa/81kDa/177kDa**

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

### Background

This gene encodes a member of the Snf2 family of helicase-like proteins. The encoded protein may play a role in DNA repair and mitochondrial function. Mutations in this gene have been associated with bone marrow failure syndrome 2. Alternatively spliced transcript variants that encode different protein isoforms have been described.

### For Research Use Only