

RBM12B Polyclonal Antibody

Catalog No. E-AB-64514

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Recombinant fusion protein of human RBM12B (NP_976324.2).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of HeLa cells using RBM12B Polyclonal Antibody at dilution of 1:3000.

Observed Mw:140kDa
Calculated Mw:118kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

RBM12B (RNA Binding Motif Protein 12B) is a Protein Coding gene. Diseases associated with RBM12B include Primary Cerebellar Degeneration and Boucher-Neuhauser Syndrome. Gene Ontology (GO) annotations related to this gene include nucleic acid binding and nucleotide binding. An important paralog of this gene is RBM12.

For Research Use Only