

Recombinant Human Carbonic Anhydrase IX/CA9 (C-Avi-6His) Biotinylated



Catalog Number:PKSH033862

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	CA9;CA-IX;Carbonic Anhydrase IX;Carbonate dehydratase IX;G250;MN;P54/58N;RCC;RCC-associated protein G250
Species	Human
Expression Host	HEK293 Cells
Sequence	Gln38-Asp414
Accession	Q16790
Calculated Molecular Weight	43.6 kDa
Observed molecular weight	60 kDa
Tag	C-Avi-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Carbonic anhydrases IX (CA IX), also known as membrane antigen MN or CA9, is a member of the carbonic anhydrase (CA) family and may be involved in cell proliferation and cellular transformation. CAs are zinc metalloenzymes that catalyze the reversible hydration of carbon dioxide ($H_2O + CO_2 = H^+ + HCO_3^-$) and thus participate in a variety of biological and physical processes. CA9 is a transmembrane enzyme expressed primarily in carcinoma cells. It is one of the best markers for hypoxia and for RCC. Appears to be a novel specific biomarker for a cervical neoplasia.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017