

Recombinant Human IL-36 Receptor Antagonist Protein/IL-36RN/IL-1F5

Catalog Number:PKSH033872

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Interleukin-36 Receptor Antagonist Protein;FIL1 Delta;IL-1-Related Protein 3;IL-1RP3;Interleukin-1 HY1;IL-1HY1;Interleukin-1 Delta;IL-1 Delta;Interleukin-1 Family Member 5;IL-1F5;Interleukin-1 Receptor Antagonist Homolog 1;IL-1ra Homolog 1;Interleukin-1-Like Protein 1;IL-1L1;IL36RN;FIL1D;IL1F5;IL1HY1;IL1L1;IL1RP3
Species	Human
Expression Host	E.coli
Sequence	Val2-Asp155
Accession	Q9UBH0
Calculated Molecular Weight	16.8 kDa
Observed molecular weight	17 kDa
Tag	None
Bioactivity	Measured by its ability to inhibit IL-36 alpha, IL-36 beta or IL-36 gamma-induced IL-8 secretion in A431 human epithelial carcinoma cells. The ED ₅₀ for this effect is 23.14 ng/ml.

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Human Interleukin-36 Receptor Antagonist (IL-36RN) is a secreted protein which belongs to the Interleukin 1 cytokine

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017

Recombinant Human IL-36 Receptor Antagonist Protein/IL-36RN/IL-1F5



Catalog Number:PKSH033872

family (IL-1 family). IL-36RN is predominantly expressed in keratinocytes but not in fibroblasts, endothelial cells or melanocytes. IL-36RN is also detected in the spleen, brain leukocyte and macrophage cell types. Increased in lesional psoriasis skin. IL-36RN is a highly and a specific antagonist of the IL-1 receptor-related protein 2-mediated response to Interleukin 1 family member 9 (IL1F9). Dysregulated expression of novel agonistic and antagonistic IL-1 family member ligands can promote cutaneous inflammation, revealing potential novel targets for the treatment of inflammatory skin disorders. Human and mouse IL-36RN share 90% sequence identity.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: www.elabscience.com

Email: techsupport@elabscience.com