

Recombinant Human GAS6 (C-6His)

Catalog Number:PKSH033928



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	AXLLG;AXLLGAXL stimulatory factor;AXSFAXL receptor tyrosine kinase ligand;Gas6;GAS-6;growth arrest-specific 6;growth arrest-specific protein 6
Species	Human
Expression Host	HEK293 Cells
Sequence	Ala31-Ala678
Accession	Q14393-2
Calculated Molecular Weight	72.7 kDa
Observed molecular weight	80-90 kDa
Tag	C-His
Bioactivity	Immobilized Human AXL-His at 10µg/ml (100 µl/well) can bind Human GAS6-His : Biotinylated by NHS-biotin prior to testing. The ED50 of Recombinant Human GAS6-His is 0.04466 ug/ml.


Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Shipping	This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.
Formulation	Supplied as a 0.2 µm filtered solution of PBS, 10% Glycerol, pH 7.4.
Reconstitution	Not Applicable

Data



> 95 % as determined by reducing SDS-PAGE.



Immobilized Human AXL-His at 10µg/ml (100 µl/well) can bind Human GAS6-His : Biotinylated by NHS-biotin prior to testing. The ED50 of Recombinant Human GAS6-His is 0.04466 ug/ml.

Background

GAS6 (Growth arrest-specific protein 6) is also known as AXL receptor tyrosine kinase ligand, AXLLG, is a multimodular protein that is up-regulated by a wide variety of cell types in response to growth arrest. Gas6 binds and induces signaling through the receptor tyrosine kinases Axl, Dtk, and Mer whose signaling is implicated in cell growth and survival, cell adhesion and cell migration. GAS6/AXL signaling plays a role in various processes such as endothelial cell survival during acidification by preventing apoptosis, optimal cytokine signaling during human natural killer cell

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017

Recombinant Human GAS6 (C-6His)

Catalog Number:PKSH033928



development, hepatic regeneration, gonadotropin-releasing hormone neuron survival and migration, platelet activation, or regulation of thrombotic responses.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: www.elabscience.com

Email: techsupport@elabscience.com