

Recombinant Human CD81 (N-FC)

Catalog Number:PKSH034009



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	CD81 antigen;CD81 molecule;CD81;CVID6;TAPA1;TAPA-1;TSPAN28
Species	Human
Expression Host	HEK293 Cells
Sequence	Phe113-Lys201
Accession	P60033
Calculated Molecular Weight	36.1 kDa
Observed molecular weight	37 kDa
Tag	N-Fc

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

CD81, also known as TAPA-1 and Tetraspanin-28, belongs to the transmembrane 4 superfamily, also known as the tetraspanin family. Members of this family mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. CD81 is a widely expressed cell-surface protein involved in an astonishing variety of biologic responses. CD81 associates with a wide range of membrane proteins including CD151, Tfr2, LDL R, PCSK9, Glypican 3, IFITM1, IGSF8/CD316, FPRP, and complexes of CD19-CD21. It is related to adhesion, morphology, activation, proliferation, and differentiation of B, T, and other cells. CD81 additionally functions as a receptor for the E2 glycoprotein of hepatitis C virus. The CD81-E2 interaction inhibits NK cell cytolytic activity, provides a co-stimulatory signal to T cells, and inhibits the maturation of plasmacytoid dendritic cells.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017