

DLGAP5 Polyclonal Antibody

Catalog No. E-AB-64727

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human DLGAP5 (NP_055565.3).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

IF 1:50-1:200

Data



Immunofluorescence analysis of HeLa cells using DLGAP5 Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

DLGAP5 (DLG Associated Protein 5) is a Protein Coding gene. Among its related pathways are Aurora A signaling. Gene Ontology (GO) annotations related to this gene include phosphoprotein phosphatase activity. An important paralog of this gene is DLGAP2.

For Research Use Only