

TWF2 Polyclonal Antibody

Catalog No. E-AB-64768

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description


Reactivity	Human, Mouse
Immunogen	Recombinant protein of human TWF2
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB 1:500 - 1:2000 IP


1:20 - 1:50

Data



Western blot analysis of extracts of various cell lines with TWF2 Polyclonal Antibody

Observed Mw:40kDa
Calculated Mw:39kDa



Immunoprecipitation analysis of 200ug extracts of 293T cells with TWF2 Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

The protein encoded by this gene was identified by its interaction with the catalytic domain of protein kinase C-zeta. The encoded protein contains an actin-binding site and an ATP-binding site. It is most closely related to twinfilin (PTK9), a conserved actin monomer-binding protein.

For Research Use Only