

URM1 Polyclonal Antibody

Catalog No. E-AB-64834

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Mouse,Rat
Immunogen	Recombinant fusion protein of human URM1 (NP_001252511.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB	1:500-1:2000
-----------	--------------

Data



Western blot analysis of extracts of Rat eye using URM1 Polyclonal Antibody at dilution of 1:3000.

Observed Mw:11kDa

Calculated Mw:7kDa/11kDa/15kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

URM1 (Ubiquitin Related Modifier 1) is a Protein Coding gene. Diseases associated with URM1 include Benign Epilepsy With Centrotemporal Spikes. Among its related pathways are Gene Expression and tRNA processing. Gene Ontology (GO) annotations related to this gene include sulfurtransferase activity.

For Research Use Only