

SPSB1 Polyclonal Antibody

Catalog No. E-AB-64852

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Recombinant fusion protein of human SPSB1 (NP_079382.2).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines using SPSB1 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:31kDa
Calculated Mw:30kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

SPSB1 (SplA/Ryanodine Receptor Domain And SOCS Box Containing 1) is a Protein Coding gene. Diseases associated with SPSB1 include Congenital Muscular Dystrophy-Dystroglycanopathy A7. Among its related pathways are Innate Immune System and Class I MHC mediated antigen processing and presentation. An important paralog of this gene is SPSB4.

For Research Use Only