

MB21D1 Polyclonal Antibody

Catalog No. E-AB-64957

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description


Reactivity	Human, Mouse, Rat
Immunogen	A synthetic peptide of human MB21D1
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB 1:500 - 1:2000

IHC 1:50 - 1:100


Data




Western blot analysis of extracts of various cell lines with MB21D1 Polyclonal Antibody

Observed Mw:59kDa

Calculated Mw:49kDa/58kDa



Immunohistochemistry of paraffin-embedded human liver injury with MB21D1 Polyclonal Antibody



Immunohistochemistry of paraffin-embedded human stomach with MB21D1 Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only