

NOL12 Polyclonal Antibody

Catalog No. E-AB-64973

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Recombinant fusion protein of human NOL12 (NP_077289.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB	1:500-1:2000
-----------	--------------

Data



Western blot analysis of extracts of various cell lines using NOL12 Polyclonal Antibody at dilution of 1:3000.

Observed Mw:36kDa
Calculated Mw:24kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

NOL12 (Nucleolar Protein 12) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include rRNA binding. An important paralog of this gene is ENSG00000100101.

For Research Use Only