

## ARFRP1 Polyclonal Antibody

**Catalog No.** E-AB-64975

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


### Description

|                     |                                                      |
|---------------------|------------------------------------------------------|
| <b>Reactivity</b>   | Human,Mouse,Rat                                      |
| <b>Immunogen</b>    | Recombinant protein of human ARFRP1                  |
| <b>Host</b>         | Rabbit                                               |
| <b>Isotype</b>      | IgG                                                  |
| <b>Purification</b> | Affinity purification                                |
| <b>Conjugation</b>  | Unconjugated                                         |
| <b>Buffer</b>       | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. |

### Applications Recommended Dilution

|           |              |
|-----------|--------------|
| <b>WB</b> | 1:500-1:2000 |
|-----------|--------------|

### Data



Western blot analysis of extracts of various cell lines  
with ARFRP1 Polyclonal Antibody

**Observed Mw:25kDa**

**Calculated Mw:13kDa/17kDa/19kDa/22kDa**

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

### Background

The protein encoded by this gene is a membrane-associated GTP-ase which localizes to the plasma membrane and is related to the ADP-ribosylation factor (ARF) and ARF-like (ARL) proteins. This gene plays a role in membrane trafficking between the trans-Golgi network and endosomes. Alternatively spliced transcript variants encoding different isoforms have been identified.

### For Research Use Only