

PARP12 Polyclonal Antibody

Catalog No. E-AB-64977

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant protein of human PARP12
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications	Recommended Dilution
--------------	----------------------

WB 1:500 - 1:2000

Data



Western blot analysis of extracts of various cell lines with PARP12 Polyclonal Antibody

Observed Mw:79kDa
Calculated Mw:79kDa

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

For Research Use Only