

TMEM25 Polyclonal Antibody

Catalog No. E-AB-64990

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human TMEM25 (NP_001137506.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

IF 1:50-1:200

Data



Immunofluorescence analysis of L929 cells using TMEM25 Polyclonal Antibody at dilution of 1:100.
Blue: DAPI for nuclear staining.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

TMEM25 (Transmembrane Protein 25) is a Protein Coding gene. Diseases associated with TMEM25 include Benign Breast Phyllodes Tumor. An important paralog of this gene is CD22.

For Research Use Only