

SLC47A1 Polyclonal Antibody

Catalog No. E-AB-65034

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Rat
Immunogen	Recombinant fusion protein of human SLC47A1 (NP_060712.2).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:500-1:1000

Data



Western blot analysis of extracts of Rat liver using SLC47A1 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:62kDa

Calculated Mw:31kDa/61kDa/63kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

This gene is located within the Smith-Magenis syndrome region on chromosome 17. It encodes a protein of unknown function.

For Research Use Only