

SLC26A2 Polyclonal Antibody

Catalog No. E-AB-65035

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Recombinant protein of human SLC26A2
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines
with SLC26A2 Polyclonal Antibody

Observed Mw:83kDa
Calculated Mw:81kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

The diastrophic dysplasia sulfate transporter is a transmembrane glycoprotein implicated in the pathogenesis of several human chondrodysplasias. It apparently is critical in cartilage for sulfation of proteoglycans and matrix organization.

For Research Use Only