

GLB1L3 Polyclonal Antibody

Catalog Number:E-AB-65161

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human GLB1L3
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Formulation	PBS with 0.01% thiomersal,50% glycerol,pH7.3.

Applications Recommended Dilution

WB	1:500-1:2000
-----------	--------------

Data



Western blot analysis of extracts of various cell lines using GLB1L3 Polyclonal Antibody at 1:1000 dilution.

Observed Mw:75kDa

Calculated Mw:35kDa/74kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Predicted to enable beta-galactosidase activity. Predicted to be involved in carbohydrate metabolic process. Predicted to be active in vacuole.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017