

KLHDC2 Polyclonal Antibody

Catalog No. E-AB-65217

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Rat
Immunogen	Recombinant fusion protein of human KLHDC2 (NP_055130.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of Rat pancreas using KLHDC2 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:46kDa

Calculated Mw:28kDa/46kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Represses CREB3-mediated transcription by interfering with CREB3-DNA binding.

For Research Use Only