

PKNOX1 Polyclonal Antibody

Catalog No. E-AB-65256

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Recombinant protein of human PKNOX1
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB 1:200 - 1:2000

Data



Western blot analysis of extracts of 293T cells with
PKNOX1 Polyclonal Antibody

Observed Mw: 48kDa

Calculated Mw: 44kDa/47kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only