

MTMR4 Polyclonal Antibody

Catalog No. E-AB-65276

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human MTMR4 (NP_004678.3).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Applications Recommended Dilution

IF 1:50-1:200

Data



Immunofluorescence analysis of L929 cells using MTMR4 Polyclonal Antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunofluorescence analysis of U-2 OS cells using MTMR4 Polyclonal Antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

MTMR4 (Myotubularin Related Protein 4) is a Protein Coding gene. Diseases associated with MTMR4 include Mulibrey Nanism and Myopathy, Centronuclear, X-Linked. Among its related pathways are Signaling by GPCR and Inositol phosphate metabolism (KEGG). Gene Ontology (GO) annotations related to this gene include phosphatase activity and protein serine/threonine phosphatase activity. An important paralog of this gene is MTMR3.

For Research Use Only