

CPEB3 Polyclonal Antibody

Catalog No. E-AB-65304

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Recombinant fusion protein of human CPEB3 (NP_055727.3).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB	1:200-1:2000
-----------	--------------

Data



Western blot analysis of extracts of various cell lines using CPEB3 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:75kDa
Calculated Mw:74kDa/76kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

CPEB3 (Cytoplasmic Polyadenylation Element Binding Protein 3) is a Protein Coding gene. Among its related pathways are Oocyte meiosis. Gene Ontology (GO) annotations related to this gene include nucleic acid binding and RNA binding. An important paralog of this gene is CPEB2.

For Research Use Only