

GIT1 Polyclonal Antibody

Catalog No. E-AB-65315

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Recombinant protein of human GIT1
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB 1:200-1:2000


Data



Western blot analysis of extracts of various cell lines with GIT1 Polyclonal Antibody

Observed Mw:100kDa

Calculated Mw:19kDa/84kDa/85kDa



Western blot analysis of extracts from normal (control) and GIT1 knockout (KO) HeLa cells, using GIT1 Polyclonal Antibody at dilution of 1:3000.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only