

SDSL Polyclonal Antibody

Catalog No. E-AB-65360

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Mouse
Immunogen	Recombinant protein of human SDSL
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications	Recommended Dilution
--------------	----------------------

WB 1:200-1:2000

Data



Western blot analysis of extracts of various cell lines with SDSL Polyclonal Antibody

Observed Mw:35kDa
Calculated Mw:34kDa

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

For Research Use Only