

KIAA0391 Polyclonal Antibody

Catalog No. E-AB-62319

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Rat
Immunogen	Recombinant fusion protein of human KIAA0391 (NP_055487.2).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines using KIAA0391 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:67kDa

Calculated Mw:24kDa/56kDa/65kDa/67kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

PRORP (Protein Only RNase P Catalytic Subunit) is a Protein Coding gene. Diseases associated with PRORP include Perrault Syndrome 1 and Syndromic X-Linked Intellectual Disability Type 10. Among its related pathways are tRNA processing and Gene Expression. An important paralog of this gene is ENSG00000258790.

For Research Use Only